



Let Our Image Enhance Your Image

Superior image quality is the best way to describe this new MedRx Video Otoscope System

This user friendly system consists of a well balanced camera and probe that require no focusing. The miniature color video camera is light weight and has a precision white balance built in. The deluxe probe consists of an aluminum body and durable stainless steel conical tip that tapers down to 3mm.

The light source utilizes a 21 volt, 150 watt halogen bulb with an approximate life of 200 hours.



MedRx INCTM

www.medrx-usa.com

1200 Starkey Rd. #105, Largo, FL 33771

Toll Free 888-392-1234

727-584-9600

Fax: 727-584-9602