

All Good Things Come

In Small Packages



MedRx AVANT™ A2D+ Diagnostic Audiometer

Powerful PC-based Dual Channel Diagnostic Audiometer

- AC, BC, Speech Audiometry and Free Field Testing
- Integrated Speech Audiometry wordlists
- Optional Fully Automated Audiometry
- Optional QuickSIN™ Speech In Noise testing
- Powerful 3rd Party Counselling Tools
- Portable, small footprint (approx. 5" x 5" x 1.25")



MedRx AVANT™ REM Speech+

Complete PC-based REM/Live Speech Mapping System

- Sure Probe™ Microphone system, FFT Dual Probe Measurements
- Features DSL, NAL-NL1 and NAL-NL2 Targets
- Red (R) and Blue (L) Illuminated Ear Indicators
- Unique 3rd Party Probe Tube monitoring system
- Powerful 3rd Party Counselling Tools
- Portable, small footprint (approx. 5" x 5" x 1.25")



MedRx AVANT™ REMsp

Smallest Single Probe Live Speech Mapping System

- FFT Monaural Probe Measurements
- Features DSL, NAL-NL1 and NAL-NL2 Targets
- Lightweight and portable, very small footprint (approx. 1" x 3")

Identical software modules/functionalities to the REM Speech+ system, yet designed as a unique Pocket Sized, Plug and Play Live Speech Mapping USB Dongle



MedRx Otowave 102-4

Handheld Portable Tympanometer

- Compatible with all MedRx software modules via an infrared link
- Fast and accurate Middle Ear Measurements
- Programmable Reflex test

MedRx AVANT AUD Audiometer

PC-based Clinical High-Frequency Audiometer

- AC, BC, Speech Audiometry and Free Field Testing
- Optional High-Frequency Audiometry up to 16,000 Hz
- Optional Fully Automated Audiometry
- Optional QuickSIN™ Speech In Noise testing
- Designated Transducer Outputs, Dual Internal Amplifiers
- Powerful 3rd Party Counselling Tools

MedRx AVANT Polar HIT+ Chamber

Unique USB Powered Directional HIT-Box

- Conducts a series of 11 selectable Automated ANSI & IEC Hearing Instrument Tests
- Provides a motor driven rotating test pad that automatically generates polar plots
- Omni and Directional microphone verification
- Acoustic isolation for accurate testing
- Detailed full page reports of conducted tests

MedRx Video Companion

Bringing Video Otoscopy To Your Computer

- Superior image quality
- Battery operated LED light source
- Focus Always Feature, Compact Color Video Camera
- Optional Trackball mouse including operator foot-pedal

MedRx Ultra Vac

"The" Hearing Instrument Restoration System

- Including:*
- Pulse / Constant Power feature
 - Durable Vacuum Wand
 - Assorted wand tips and Multi-Piece Toolkit
 - Vacuum Drying Chamber
 - Vacuum Gauge and Particle Filter
 - Selectable wand (vacuum or pressure)



MedRx®

Precision Acoustic Instruments

* Most MedRx products are USB powered
** MedRx products are NOAH™ compatible
*** Lifetime Free software upgrades where applicable

www.medrx-usa.com or Call Toll Free: 888-392-1234