


SLIGHTLY LARGER
Approx.
5" x 6 1/2"
GREATLY IMPROVED

Our A2D+ Audiometer Had More To Give

Now 1 1/2" Longer To Accommodate Two Sets
Of Headphones At The Same Time



Now
With Automatic
Driver
Installation

MedRx AVANT A2D+

- Dual Channel Audiometer
- Air, Bone & Speech Functions
- Powerful 3rd Party Counseling Tools
- Optional Integrated Quick SIN™
- Optional Automated Audiometry
- Compact Design, Portable
- Small Footprint
- PC Based via USB Connection
- Runs within NOAH™ or Stand Alone

At MedRx, continuous improvement is a way of life. We listen to our customers!

Now, the AVANT™ A2D+ is truly plug and play. No need to load separate drivers during the initial device installation. Simply install the operating software, plug the device into the USB 2.0 port of your computer and you are ready to test. No longer USB port specific! The computer will recognize the audiometer immediately, making the installation process even faster and easier.

The new design also provides dual air conduction ports, allowing two separate headsets to be plugged in simultaneously. No more inconvenient plugging and unplugging of headsets.

MedRx[®]
www.medrx-usa.com

*Experience the NEW AVANT A2D+ today and
see why it's changing diagnostic audiometry worldwide*

Call 888-392-1234 or 727-584-9600