

Precision Audiometric Instruments

MedRx has created a remarkable new generation of discreet yet powerful PC-based instrumentation for Audiometry, Real Ear Measurement, Live Speech Mapping, Hearing Instrument Testing and Video Otoscopy



MedRx AVANT™ A2D+ A Powerful PC-based Dual Channel Diagnostic Audiometer

- AC, BC, Speech Audiometry & Free Field Testing
- Integrated Speech Audiometry Wordlists
- Optional Fully Automated Audiometry
- Optional QuickSIN™ Speech In Noise Testing
- Powerful 3rd Party Counseling Tools
- Portable, Small Footprint (Approximately 5" x 5" x 1.25")

MedRx AVANT™ REM Speech+ A Complete PC-based REM/Live Speech Mapping System

- Sure Probe™ Microphone System, FFT Dual Probe Measurements
- Features DSL, NAL-NL1 And NAL-NL2 Targets
- Red (R) And Blue (L) Illuminated Ear Indicators
- Unique 3rd Party Probe Tube Monitoring System
- Powerful 3rd Party Counseling Tools
- Portable, Small Footprint (Approximately 5" x 5" x 1.25")

MedRx AVANT™ REMsp Our Single Probe Live Speech Mapping System

- FFT Monaural Probe Measurements
- Features DSL, NAL-NL1 And NAL-NL2 Targets
- Lightweight And Portable, Very Small Footprint (Approximately 1" x 3")
- Input Levels Can Be Viewed In Either LTASS Or Peak Hold Mode
- Hearing Loss Simulator For 3rd Party Demonstrations

Identical software modules/functionality to the REM Speech+ system, yet designed as a unique pocket sized, plug and play Live Speech Mapping USB dongle

MedRx Otowave A Handheld Portable Tympanometer

- Compatible With All MedRx Software Modules Via An Infrared Link
- Fast And Accurate Middle Ear Measurements
- Programmable Reflex Test

MedRx AVANT "Stealth" Audiometer A PC-Based Clinical High-Frequency Audiometer

- AC, BC, Speech Audiometry And Free Field Testing
- Optional High-Frequency Audiometry Up To 16,000 Hz
- Optional Fully Automated Audiometry
- Optional QuickSIN™ Speech In Noise testing
- Designated Transducer Outputs, Dual Internal Amplifiers
- Powerful 3rd Party Counseling Tools

MedRx AVANT Polar HIT+ Chamber A Unique USB Powered Directional HIT-Box

- Conducts A Series Of 11 selectable Automated ANSI & IEC Hearing Instrument Tests
- Provides A Motor Driven Rotating Test Pad That Automatically Generates Polar Plots
- Omni And Directional Microphone Verification
- Acoustic Isolation For Accurate Testing
- Detailed Full Page Reports Of Conducted Tests

MedRx Video Companion Bringing Video Otoscopy To Your Computer

- Superior Image Quality
- Battery Operated LED Light Source
- Focus Always Feature, Compact Color Video Camera
- Optional Trackball Mouse Including Operator Foot-Pedal

MedRx Ultra Vac "The" Hearing Instrument Restoration System, Includes:

- Pulse / Constant Power Feature
- Durable Vacuum Wand
- Vacuum Drying Chamber
- Vacuum Gauge And Particle Filter
- Selectable Wand (Vacuum Or Pressure)

MedRx®

Toll Free: 888-392-1234 • www.medrx-usa.com