

Audiological Diagnostic & Video Imaging Equipment

Powerful PC-Based Instrumentation For Audiometry, Real Ear Measurement, Live Speech Mapping, Hearing Instrument Testing/Evaluation & Video Otoscopy



AVANT™ A2D+ Audiometer

Powerful PC-Based Dual Channel Diagnostic Audiometer

- AC, BC, Speech, Masking & Free Field Testing
- Integrated Speech Audiometry Wordlists
- Optional Fully Automated Audiometry
- Optional QuickSIN™ Speech In Noise Testing
- Powerful 3rd Party Counseling Tools
- Small Footprint, Approx 5" x 5" x 1.25"



AVANT REM Speech+

Complete Binaural Live Speech Mapping, Real Ear Measurement System

- Dual Sure Probe™ Microphone System
- USB Powered
- Results Can Be Viewed In Either LTASS Or Peak Hold Mode
- Features DSL (I/O), NAL-NL1 And NAL-NL2 Targets
- Hearing Loss Simulator For 3rd Party Demonstrations
- Small Footprint, Approx 5" x 5" x 1.25"



AVANT "Stealth" Audiometer

PC-Based Clinical Audiometer

- AC, BC, Speech, Masking And Free Field Testing
- Designated Transducer Outputs
- Dual Internal Amplifiers
- Powerful 3rd Party Counseling Tools
- Optional: Fully Automated Audiometry, QuickSIN™ Speech In Noise Testing, High-Frequency Audiometry

AVANT REMsp

Monaural Live Speech Mapping, Real Ear Measurement System

- USB Powered
- Results Can Be Viewed In Either LTASS Or Peak Hold Mode
- Features DSL (I/O), NAL-NL1 And NAL-NL2 Targets
- Hearing Loss Simulator
- Approx 1" x 3" x 0.5"



AVANT Polar HIT+ Chamber

Unique USB Powered Directional Hearing Instrument Test Box

- Conducts 11 Selectable Automated ANSI & IEC Hearing Instrument Tests
- Acoustic Isolation For Accurate Testing
- Provides A Motor Driven Rotating Test Pad To Automatically Generate Polar Plots
- Omni And Directional Microphone Verification

Otowave Tymp

Handheld Portable Tympanometer

- Compatible With All MedRx Software Modules Via An Infrared Link
- Programmable Reflex Testing
- Fast And Accurate Middle Ear Measurement



Video Companion

Bringing Video Otoscopy To Your Computer

- Superior Image Quality
- Battery Operated LED Or Halogen Light Source
- Compact, Color Video Camera
- Focus Always Feature
- Optional Trackball Mouse Including Operator Foot-Pedal

Ultra Vac

Hearing Instrument Restoration System, Includes:

- Pulse / Constant Power Feature
- Durable Vacuum Wand
- Vacuum Drying Chamber
- Selectable Wand (Vacuum Or Pressure)

MedRx Products Are NOAH™ And TIMS® Compatible

OT-S-ADH169-1

MedRx®

www.medrx-usa.com • Toll Free 888-392-1234