



Diagnostic & Testing Instrumentation

COMPUTER REQUIREMENTS

MedRx **Minimum** Computer Specs:

- Windows®-PC compatible computer
- Intel™ i5 Dual Core, 2.0 GHz or better
- 4 GB RAM
- 20 GB free hard drive space
- Available 2.0 USB port
- Windows 7, 8 or 10 Professional (32 or 64-bit)

MedRx **Recommended** Computer Specs:

- Windows®-PC Computer
- Intel™ i5 Dual Core, 3.2 GHz or better
- 8 GB RAM or more
- Available USB 2.0
- Graphics Adapter with 2GB Dedicated Video Memory
- 50 GB or more free hard drive space
- DVD-ROM Drive
- High Speed Internet Connection
- Windows 10 Professional 64-bit

MedRx®